

**Agilent Human 14  
Multiple Affinity  
Removal System  
Columns for the  
Fractionation of  
High-Abundant  
Proteins from Human  
Proteomic Samples**

Instructions



Agilent Technologies

# General Information

## Introduction


The Agilent Human 14 Multiple Affinity Removal System consists of a Multiple Affinity Removal System Column based on affinity interactions and optimized buffers for sample loading, washing, eluting, and regenerating. It is specifically designed to fractionate 14 high-abundant proteins from human biological fluids such as plasma, serum, and cerebral spinal fluid (CSF). This technology enables removal of albumin, IgG, antitrypsin, IgA, transferrin, haptoglobin, fibrinogen, alpha2-macroglobulin, alpha1-acid glycoprotein, IgM, apolipoprotein AI, apolipoprotein AII, complement C3, and transthyretin in a single step. The targeted high-abundant proteins are simultaneously removed when crude biological samples are passed through the column. Selective immunodepletion provides an enriched pool of low-abundant proteins for downstream proteomics analysis.

Specific removal of 14 high-abundant proteins depletes approximately 94% of total protein mass from human serum. The low-abundant proteins in the flow-through fractions can be studied. Removal of high-abundant proteins enables improved resolution and dynamic range for one-dimensional gel electrophoresis (1DGE), two-dimensional gel electrophoresis (2DGE), and liquid chromatography/

mass spectrometry (LC/MS). The collected flow-through fractions may need to be concentrated dependent upon the downstream applications.

# Human 14 Multiple Affinity Removal System

The Agilent Human 14 Multiple Affinity Removal System is specially designed to fractionate the high-abundant proteins from the low-abundant proteins, allowing improved analysis of both fractions. (see Figure 1).



**Figure 1** The Multiple Affinity Removal System

# Product Description

**Table 1 Human 14 Multiple Affinity Removal System Column and Starter Reagent Kit**

<b>Product no.</b>	<b>Product name</b>	<b>Product description</b>
5188-6557	4.6 × 50 mm, affinity column	Remove human albumin, IgG, antitrypsin, IgA,
5188-6558	4.6 × 100 mm, affinity column	transferrin, haptoglobin and fibrinogen,
5188-6559	10 × 100 mm, affinity column	alpha2-macroglobulin, alpha1-acid glycoprotein, IgM, apolipoprotein AI, apolipoprotein AII, complement C3, transthyretin
5185-5987	Buffer A, 1 L	Ready-to-use, optimized buffer for loading, washing, and equilibrating column
5185-5988	Buffer B, 1 L	Ready-to-use, optimized buffer for elution of bound proteins from column
5185-5990	Spin filters, 0.22 µm, 1 pack of 25	For sample cleanup before loading column
5185-5991	Concentrators, 5 kDa, MWCO, 1 pack of 25	For concentrating flow-through fractions
5185-5989	Human serum albumin	Dilute standard for checking column capacity

Product no.	Product name	Product description
5185-5995	Replacement frit for 4.6 mm ID column (1 each)	Replaced when column inlet frit is clogged.
5021-8814	Replacement frit for 10 mm ID column (1 each)	Replace when column inlet frit is clogged.
5185-5986	<b>Starter Reagent Kit*</b> Buffer A: 2 × 1 L Buffer B: 1 L Spin filters 0.22 µm: 2 packs of 25 Protein concentrators: 1 pack of 25	

\* Under normal use conditions, the kit should last for approximately 200 injections for 4.6 × 50 mm columns, and 100 injections for 4.6 × 100 mm columns.

**CAUTION**

Do not expose columns to organic solvents (like alcohols, acetonitrile, etc.), strong oxidizers, acids, or reducing agents, and other protein-denaturing agents. Before attaching the column, purge the LC system and run two method blank injections according to protocol to ensure all lines and sample loops are free of organic solvents.

**NOTE**

For LC systems shared with other chemical applications, be sure to purge the LC system first with isopropyl alcohol, and then extensively with water (approximately 1 hour). After purging, proceed with protocol.

# Full Protocol for 4.6 x 50 mm Column

(Column capacity: 20  $\mu$ L human serum/plasma\*)

- 1 Set up Buffer A and Buffer B as the only mobile phases.
- 2 Purge lines with Buffer A and Buffer B at a flow rate of 1.0 mL/min for 10 min.
- 3 Set up LC timetable (see Table 1 for details) and run two method blanks by injecting 100  $\mu$ L of Buffer A **without a column**.

## NOTE

Ensure proper sample loop size in autosamplers.

- \* Consult column certificate of analysis to verify your column capacity.

It is the user's responsibility to adjust loading volume to compensate for the fact that the concentrations of some high-abundant proteins can vary widely depending on the sample origin (source). Several plasma proteins such as alpha1-antitrypsin, haptoglobin, fibrinogen, IgG, and alpha1-acid glycoprotein rise several folds in response to stress, infection, inflammation, or tissue necrosis and are known as acute phase reactants (Henry, J.B. 1996 - Clinical Diagnosis and Management by Laboratory Methods). Capacity testing was done with pooled normal plasma from Rockland Immunochemicals (D519-04).



**Table 2 LC method for 4.6 x 50 mm column**

Solvent A: Buffer A

Pressure limits: 60 bar


Solvent B: Buffer B

**LC Timetable**

	<b>Time</b>	<b>%B</b>	<b>Flow rate</b>	<b>Max. pressure</b>
1	0.00	0.00	0.125	60
2	9.50	0.00	0.125	60
3	9.51	0.00	1.000	60
4	11.50	0.00	1.000	60
5	11.51	100.00	1.000	60
6	16.00	100.00	1.000	60
7	16.01	0.00	1.000	60
8	25.00	0.00	1.000	60

- 4** Attach column and equilibrate it in Buffer A for 4 min at a flow rate of 1 mL/min at room temperature.

- 5 Dilute human serum/plasma\* 4 times (for example: 20  $\mu$ L human plasma with 60  $\mu$ L of Buffer A†).
- 6 Remove particulates with a 0.22  $\mu$ m spin filter; 1 min at 16,000 g.
- 7 Inject 80  $\mu$ L of the diluted serum/plasma at a flow rate of 0.125 mL/min (injected volume should be adjusted to your column capacity).
- 8 Collect flow-through fraction (like that which appears between 5.0–7.5 min in Figure 2) and store collected fractions at -20 °C if not analyzed immediately.



**Figure 2** Representative chromatogram of 4.6 x 50 mm column.

- \* Protocol may be applied to other human biological fluids like CSF.
- † Addition of protease inhibitors in Buffer A for sample dilution helps prevent protein degradation.

- 9** Elute bound proteins from the column with Buffer B (elution buffer) at a flow rate of 1 mL/min for 5.5 min.
- 10** Regenerate column by equilibrating it with Buffer A for an additional 9.0 min at a flow rate of 1 mL/min.
- 11** Store column after equilibrating with Buffer A at 2–8 °C (35–46 °F) in a refrigerator. DO NOT FREEZE THE COLUMN.
- 12 Analyze.** Analyze flow-through containing the low-abundant proteins. For 1D-SDS-PAGE, an aliquot of the flow-through fraction may be used directly. For IEF, 2DGE, and MS analysis of the flow-through fraction, it is necessary to buffer exchange/desalt to an appropriate buffer. The 5KDa MWCO spin concentrators (part number 5185-5991) may be used for buffer exchange and concentration. Alternatively, for automated desalting and concentration, the Agilent mRP-C18 column (part number 5188-5231) may be used according to published methods (Agilent Technologies, publication 5989-2506EN).

# Full Protocol for 4.6 x 100 mm Column

(Column capacity: 40  $\mu$ L human serum/plasma\*)

- 1 Set up Buffer A and Buffer B as the only mobile phases.
- 2 Purge LC system with Buffer A and Buffer B at a flow rate of 1.0 mL/min for 10 min **without a column**.
- 3 Set up LC timetable (see Table 3 for details) and run two method blanks by injecting 200  $\mu$ L of Buffer A **without a column**.

## NOTE

Ensure proper sample loop size in autosamplers.

- \* Consult column certificate of analysis to verify your column capacity.

It is the user's responsibility to adjust loading volume to compensate for the fact that the concentrations of some high-abundant proteins can vary widely depending on the sample origin (source). Several plasma proteins such as alpha1-antitrypsin, haptoglobin, fibrinogen, IgG, and alpha1-acid glycoprotein rise several folds in response to stress, infection, inflammation, or tissue necrosis and are known as acute phase reactants (Henry, J.B. 1996 - Clinical Diagnosis and Management by Laboratory Methods). Capacity testing was done with pooled normal plasma from Rockland Immunochemicals (D519-04).

**Table 3 LC method for 4.6 x 100 mm column**

Solvent A: Buffer A

Pressure limits: 60 bar


Solvent B: Buffer B

**LC Timetable**

	<b>Time</b>	<b>%B</b>	<b>Flow rate</b>	<b>Max. pressure</b>
1	0.00	0.00	0.125	60
2	18.00	0.00	0.125	60
3	18.01	0.00	1.000	60
4	20.00	0.00	1.000	60
5	20.01	100.00	1.000	60
6	27.00	100.00	1.000	60
7	27.01	0.00	1.000	60
8	38.00	0.00	1.000	60

- 4 Attach column and equilibrate it in Buffer A for 4 min at a flow rate of 1 mL/min at room temperature.

- 5 Dilute human serum/plasma\* 4 times (for example: 40  $\mu$ L human plasma with 120  $\mu$ L of Buffer A†).
- 6 Remove particulates with a 0.22  $\mu$ m spin filter; 1 min at 16,000 g.
- 7 Inject 160  $\mu$ L of the diluted serum/plasma at a flow rate of 0.125 mL/min (injected volume should be adjusted to your column capacity).
- 8 Collect flow-through fraction (like that which appears between 11–15 min in Figure 3) and store collected fractions at -20 °C if not analyzed immediately.



**Figure 3** Representative chromatogram of 4.6 x 100 mm column.

- \* Protocol may be applied to other human biological fluids like CSF.
- † Addition of protease inhibitors in Buffer A for sample dilution helps prevent protein degradation.

- 9** Elute bound proteins from the column with Buffer B (elution buffer) at a flow rate of 1 mL/min for 7.0 min.
- 10** Regenerate column by equilibrating it with Buffer A for 11.0 min at a flow rate of 1 mL/min.
- 11** Store column after equilibrating with Buffer A at 2–8 °C (35–46 °F) in a refrigerator. **DO NOT FREEZE THE COLUMN.**
- 12 Analyze.** Analyze flow-through containing the low-abundant proteins. For 1D-SDS-PAGE, an aliquot of the flow-through fraction may be used directly. For IEF, 2DGE, and MS analysis of the flow-through fraction, it is necessary to buffer exchange/desalt to an appropriate buffer. The 5KDa MWCO spin concentrators (part number 5185-5991) may be used for buffer exchange and concentration. Alternatively, for automated desalting and concentration, the Agilent mRP-C18 column (part number 5188-5231) may be used according to published methods (Agilent Technologies, publication 5989-2506EN).

# Full Protocol for 10 x 100 mm Column

(Column capacity: 250 µL human serum/plasma\*)

- 1 Set up Buffer A and Buffer B as the only mobile phases.
- 2 Purge LC system with Buffer A and Buffer B at a flow rate of 1.0 mL/min for 10 min **without a column**.
- 3 Set up LC timetable (see Table 2 for details) and run two method blanks by injecting 200 µL of Buffer A **without a column**.

## NOTE

Ensure proper sample loop size in autosamplers.

- \* Consult column certificate of analysis to verify your column capacity.

It is the user's responsibility to adjust loading volume to compensate for the fact that the concentrations of some high-abundant proteins can vary widely depending on the sample origin (source). Several plasma proteins such as alpha1-antitrypsin, haptoglobin, fibrinogen, IgG, and alpha1-acid glycoprotein rise several folds in response to stress, infection, inflammation, or tissue necrosis and are known as acute phase reactants (Henry, J.B. 1996 - Clinical Diagnosis and Management by Laboratory Methods). Capacity testing was done with pooled normal plasma from Rockland Immunochemicals (D519-04).



**Table 4 LC Method for 10 x 100 mm Column**

Solvent A: Buffer A

Pressure limits: 160 bar


Solvent B: Buffer B

**LC Timetable**

	<b>Time</b>	<b>%B</b>	<b>Flow rate</b>	<b>Max. pressure</b>
1	0.00	0.00	0.50	160
2	20.00	0.00	0.50	160
3	20.01	0.00	1.000	160
4	22.50	0.00	1.000	160
5	22.51	100.00	3.000	160
6	30.00	100.00	3.000	160
7	30.01	0.00	3.000	160
8	40.00	0.00	3.000	160

- 4** Attach column and equilibrate it in Buffer A for 4 min at a flow rate of 1 mL/min at room temperature.

- 5 Dilute human serum/plasma\* 4 times (for example: 200  $\mu$ L human plasma with 600  $\mu$ L of Buffer A<sup>†</sup>).
- 6 Remove particulates with a 0.22  $\mu$ m spin filter; 1 min at 16,000 g.
- 7 Inject 800  $\mu$ L of the diluted plasma at a flow rate of 0.5 mL/min (injected volume should be adjusted to your column capacity).
- 8 Collect flow-through fraction (like that which appears between 10.5–14.5 min in Figure 4) and store collected fractions at -20 °C if not analyzed immediately.



**Figure 4** Representative chromatogram of 4.6 x 100 mm column.

- \* Protocol may be applied to other human biological fluids like CSF.
- † Addition of protease inhibitors in Buffer A for sample dilution helps prevent protein degradation.

- 9** Elute bound proteins from the column with Buffer B (elution buffer) at a flow rate of 3 mL/min for 7.5 min.
- 10** Regenerate column by equilibrating it with Buffer A for 10.0 min at a flow rate of 3 mL/min.
- 11** Store column after equilibrating with Buffer A at 2–8 °C (35–46 °F) in a refrigerator. DO NOT FREEZE THE COLUMN.
- 12 Analyze.** Analyze flow-through containing the low-abundant proteins. For 1D-SDS-PAGE, an aliquot of the flow-through fraction may be used directly. For IEF, 2DGE, and MS analysis of the flow-through fraction, it is necessary to buffer exchange/desalt to an appropriate buffer. The 5KDa MWCO spin concentrators (part number 5185-5991) may be used for buffer exchange and concentration. Alternatively, for automated desalting and concentration, the Agilent mRP-C18 column (part number 5188-5231) may be used according to published methods (Agilent Technologies, publication 5989-2506EN).

# Troubleshooting

It is the user's responsibility to adjust loading volume to compensate for the fact that the concentrations of some high-abundant proteins can vary widely depending on the sample origin (source). Several plasma proteins such as alpha1-antitrypsin, haptoglobin, fibrinogen, IgG, and alpha1-acid glycoprotein rise several folds in response to stress, infection, inflammation, or tissue necrosis and are known as acute phase reactants (Henry, J.B. 1996 - Clinical Diagnosis and Management by Laboratory Methods).

Review the following information for troubleshooting your experiments.

## **High backpressure**

A clogged inlet frit may result in high backpressure, distorted peak shape, and diminished column lifetime. To prevent these problems, remove particulates from samples with a spin filter before loading. If clogging occurs, replace plugged inlet frit (Part numbers: 5185-5995 for 4.6 mm ID columns and 5021-8814 for 10 mm ID columns).

## **No bound fraction peak**

Bound proteins can only be removed from the column by eluting with Buffer B. Check LC timetable to ensure enough column time exposure to Buffer B, for complete removal of bound proteins.

## **Abnormal peak height**

Approximately 94% of serum/plasma proteins will be removed as the bound fraction. The peak height of the bound fraction is expected to be much greater than that of the flow-through fraction. If this order is reversed, two possibilities may be checked:

- a** Column may not have been regenerated well enough from previous runs, resulting in lost capacity. To correct this, elute bound proteins with Buffer B for 3 additional minutes and re-equilibrate the column with Buffer A.
- b** Check for signs of biological growth in the Buffer A reservoir. Replace with fresh Buffer A for optimal column performance.

## **Breakthrough of high-abundant proteins in flow-through fraction**

The majority (95–99%) of targeted proteins should be depleted from the sample. If the flow-through fraction contains unusually high levels of the 14 depleted proteins, two possibilities may be checked:

- a** The column capacity may have been exceeded. Reduce the sample load per injection.
- b** Serum/plasma levels of the targeted proteins may be unusually high.

# Column Specifications

**Table 5 Column Descriptions**

<b>Parameter</b>	<b>Description</b>
Part number 5188-6557	4.6 mm × 50 mm (0.83 mL)
Part number 5188-6558	4.6 mm × 100 mm (1.66 mL)
Column body material	Polyetheretherketone (PEEK)
Part Number 5188-6559	10 mm × 100 mm (7.85 mL)
Column body material	PEEK encapsulated in stainless steel
End-fitting materials	PEEK with 2 µm frits
Column capacity*	0.83 mL column: 20 µL human plasma 1.66 mL column: 40 µL human plasma 7.85 mL column: 250 µL human plasma
Max. pressure	60 or 160 bar

**Table 5 Column Descriptions**

<b>Parameter</b>	<b>Description</b>
Operating temperature	18–25 °C
Column packing material	Affinity resin
Immobilized ligands	Affinity ligands to human albumin, IgG, antitrypsin, IgA, transferrin, haptoglobin and fibrinogen, alpha2-macroglobulin, alpha1-acid glycoprotein, IgM, apolipoprotein AI, apolipoprotein AII, complement C3, transthyretin
Flow rate range	0.125–3.0 mL/min
Shipping solution	Buffer A with < 0.02% sodium azide
Shipping temperature	2–8 °C (35–46 °F)
Storage temperature	2–8 °C (35–46 °F)

\* For exact column capacity, consult your column certificate of analysis.

## Safety Issues

When preparing biological samples using any Agilent Multiple Affinity Removal System Column, follow general guidelines for handling biological materials and wear protective eyewear and gloves.

## Recommendations

- **Sample dilution**  
It is not recommended to load crude serum/plasma directly onto the column. Follow instructions for serum/plasma dilution (four times dilution with Buffer A). Addition of protease inhibitors in Buffer A for sample dilution helps prevent protein degradation.
- **Sample cleanup**  
Human serum/plasma may contain particulate materials that can be removed by a quick spin using a 0.22- $\mu\text{m}$  spin filter.
- **Column performance**  
Always use Buffers A and B as the mobile phases according to the LC protocols. Do not expose the columns to solvents other than Buffers A and B.
- **Column storage**  
After equilibrating Agilent Multiple Affinity Removal System Columns with Buffer A, always store them in a refrigerator at 2–8 °C (35–46 °F) when not in use to minimize loss in column capacity.



- **Lyophilization of flow-through fractions**  
Buffer exchange to a volatile buffer (for example, ammonium bicarbonate) is recommended prior to lyophilization due to high salt concentration in Buffer A.
- **Multiple Affinity Removal System**  
The Agilent Multiple Affinity Removal System is specially designed to allow for close study of low-abundant proteins present in flow-through fractions (see Figure 1).

**CAUTION**

Do not expose columns to organic solvents (like alcohols, acetonitrile, etc.), strong oxidizers, acids, or reducing agents, and other protein denaturing agents.

**WARNING**

**For RESEARCH USE ONLY. This product is NOT TO BE USED AS AN *IN-VITRO* DIAGNOSTIC.**

**Storage:** Upon its receipt and when you are not using it, store the column with the end-caps tightly sealed in a refrigerator at 2–8 °C (35–46 °F). **DO NOT FREEZE THE COLUMN.**

For more information or technical assistance, please call toll free: 1-800-227-9770 or visit our Web site at: [www.agilent.com/chem/bioreagents](http://www.agilent.com/chem/bioreagents).

## Related Agilent Products

Other related Agilent products include the following:

- 5188-6560 Agilent Multiple Affinity Removal System Spin Cartridge Hu-14**, spin cartridge that depletes 14 high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, antitrypsin, fibrinogen, alpha2-macroglobulin, alpha1-acid glycoprotein, IgM, apolipoprotein AI, apolipoprotein AII, complement C3, and transthyretin) from human serum/plasma samples, 8–10  $\mu$ L serum/plasma capacity per use
- 5188-6409 Agilent High Capacity Multiple Affinity Removal System Column Hu-7HC, 4.6 X 50 mm**, LC column that depletes seven high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, antitrypsin, and fibrinogen) from human serum or plasma samples, 30–35  $\mu$ L serum/plasma capacity per injection
- 5188-6410 Agilent High Capacity Multiple Affinity Removal System Column Hu-7HC, 4.6 X 100 mm**, LC column that depletes seven high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, antitrypsin, and fibrinogen) from human serum or plasma samples, 60–70  $\mu$ L serum/plasma capacity per injection

- 5188-6411 Agilent High Capacity Multiple Affinity Removal System Column Hu-7HC, 10 X 100 mm**, LC column that depletes seven high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, antitrypsin, and fibrinogen) from human serum or plasma samples, 250–300  $\mu$ L serum/plasma capacity per injection
- 5188-6408 Agilent High Capacity Multiple Affinity Removal System Spin Cartridge Hu-7HC**, spin cartridge that depletes seven high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, antitrypsin, and fibrinogen) from human serum or plasma samples, 12–14  $\mu$ L serum/plasma capacity per use
- 5188-5332 Agilent High Capacity Multiple Affinity Removal System Column Hu-6HC, 4.6 X 50 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 30–40  $\mu$ L serum capacity per injection
- 5188-5333 Agilent High Capacity Multiple Affinity Removal System Column Hu-6HC, 4.6 X 100 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 60–80  $\mu$ L serum capacity per injection

- 5188-5336 Agilent High Capacity Multiple Affinity Removal System Column Hu-6HC, 10 X 100 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 300–325  $\mu$ L serum capacity per injection
- 5188-5341 Agilent High Capacity Multiple Affinity Removal System Spin Cartridge Hu-6HC**, spin cartridge that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 14–16  $\mu$ L serum capacity per use
- 5185-5984 Multiple Affinity Removal System Column Hu-6, 4.6  $\times$  50 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 15–20  $\mu$ L serum capacity per injection
- 5185-5985 Multiple Affinity Removal System Column Hu-6, 4.6  $\times$  100 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 30–40  $\mu$ L serum capacity per injection

- 5188-2714 Multiple Affinity Removal System Column Hu-6, 10 × 100 mm**, LC column that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 140–190  $\mu$ L serum capacity per injection
- 5188-5230 Multiple Affinity Removal System Spin Cartridge Hu-6**, spin cartridge that depletes six high-abundant proteins (albumin, IgG, IgA, transferrin, haptoglobin, and antitrypsin) from human serum/plasma samples, 7–10  $\mu$ L serum capacity per use
- 5188-5217 Multiple Affinity Removal System Column Ms-3, 4.6 × 50 mm**, LC column that depletes three high-abundant proteins (albumin, IgG, and transferrin) from mouse serum/plasma samples, 37–50  $\mu$ L serum capacity per injection
- 5188-5218 Multiple Affinity Removal System Column Ms-3, 4.6 × 100 mm**, LC column that depletes three high-abundant proteins (albumin, IgG, and transferrin) from mouse serum/plasma, 75–100  $\mu$ L serum capacity per injection
- 5188-5289 Multiple Affinity Removal System Spin Cartridge Ms-3**, spin cartridge that depletes three high-abundant proteins (albumin, IgG, and transferrin) from mouse serum/plasma samples, 25–30  $\mu$ L serum capacity per use

- 5185-5986 Multiple Affinity Removal System Reagent Kit**, starter reagent kit containing buffers, spin filters, and spin concentrators for use with Multiple Affinity Removal System LC Columns
- 5188-5254 Starter Reagent Kit for Spin Cartridges\***,  
Buffer A: 1 L  
Buffer B: 1 L  
Spin filters 0.22  $\mu\text{m}$ : 2 packs of 25  
Protein concentrators: 1 pack of 25  
Luer-Lock adapters: 1 pack of 2  
5-mL plastic Luer-Lock syringes: 1 pack of 2  
1.5-mL microtubes: 6 packs of 100  
Spin cartridge extra caps and plugs, 1 pack of 6 each  
Teflon Luer-Lock needles, 1 pack of 10
- 5188-5231 mRP-C18 High Recovery Protein Fractionation and Desalting Column**, see [www.agilent.com/chem/bioreagents](http://www.agilent.com/chem/bioreagents) for more details

\* Under normal conditions, the kit should last for approximately 200 spin cartridge uses.

For more information or technical assistance, please call toll free: 1-800-227-9770 or visit our Web site at: [www.agilent.com/chem](http://www.agilent.com/chem).



©Agilent Technologies, Inc. 2007

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

Printed in the USA  
May 2007  
5969-1596



**Agilent Technologies**